







Order form for custom-made Compression stump shrinkers



Instructions for precise measuring of Juzo® compression stump shrinkers

Please note:

To ensure a correct pressure distribution, please measure the first circumference 5 cm proximal to the end of the stump. Record this in the relevant place on the measurement sheet. As per the schematic, you need to provide **all** measurements above this according to the style you have selected.

Please ensure that you complete all the necessary measurements when a hip attachment or body part is included.

The length measurements are to be taken from the end of the stump.

Measure the lengths in a straight line, do not put the measuring tape along the stump.


If necessary, use a measuring board to obtain accurate length measurements.

In the case of **compression stump shrinker pantyhose**, the following options are available:

- One stocking leg for stump treatment: Please send both measuring sheets “stump shrinker” and “stockings/pantyhose flat knitted” and include a suitable mark.
- Both stocking legs for stump treatment. Determine all the measurements up to “T”.
- One stocking leg for stump treatment, the second with a leg extension. In this case please specify also the circumference and length measurements for “F” and “G”.

Compression stump shrinker with body part (KT Juzo® Expert body part – height 23 cm)

Available sizes	1	2	3	4	5
To suit waist	90	100	110	120	135



For the selection of the appropriate garment, we recommend the brochure Juzo® Expert Special-purpose medical garments.

Quantity: piece

Former order no./dd:

Please check:

Compression classes

1

2

circular knitted (seamless)

Juzo® Dynamic¹⁾

Sesame Almond Black pepper

Juzo® Dynamic Silver – colour: Almond

¹⁾ If nothing is noted, we deliver colour Almond

flat knitted (with seam)

Juzo® Expert²⁾

White Light beige Beige Dark brown
 Red Pink Green Blue
 Dark blue Grey Black

Juzo® Expert Silver – colour: Beige

Juzo® Expert Cotton – colour: Beige

²⁾ If nothing is noted, we deliver colour Beige

3511

3512

3511

3512

3021

3022

3021

3022

3021

3022

Compression stump shrinker

with silicone border

with hip attachment (please mention circumference "cT")

stump: left right both-sided, to be worn as a pair

with body part Juzo® Expert (open crotch)

stump: left right

body part custom-made

standard body part – colour: Beige

size: length: normal short

Compression Stump shrinker "Pirogoff"

YD after Pirogoff-amputation

Compression stump shrinker pantyhose

stump: left right both-sided

open crotch with fly

Treatment of the 2nd leg:

compression stocking from "cA" (please enclose measuring sheet)

leg extension with top welt at "cF"

Circumferences in cm


Length measurements in cm

patient's left

patient's right

patient's left

patient's right



Measure circumferences beginning 5 cm from the distal end of the stump.

lST
lSH
lSK	}	}
lSG		
lSF
lSE
lSD
lSC
lSB ¹
lSB
lSY

S = distal end of the stump, measure all lengths from "S"

From:

Contact person:

Special requests:

.....

.....

.....

.....

.....

.....

Length measurements in cm

Circumferences in cm

		Standard sizes					
		I XS	II S	III M	IV L	V XL	VI XXL
<p>Compression stump shrinkers above knee Art. no. 3511, 3512</p>	cG	42-45	45-49	49-53	53-57	57-61	61-65
	cF	30-32	32-35	35-38	38-41	41-44	44-47
<p>Compression stump shrinkers below knee Art. no. 3511, 3512</p>	cF	40-42	42-45	45-48	48-51	51-54	54-57
	cE	30-32	32-35	35-38	38-41	41-44	44-47
	cC	28-30	30-33	33-36	36-39	39-42	42-45

Julius Zorn GmbH
 Juliusplatz 1
 86551 Aichach
 Germany

Telephone 00 49/82 51/90 13 00
 Telefax 00 49/82 51/90 13 21
 E-mail: order@juzo.de

www.juzo.com

Juzo UK Ltd
 South Court
 Sharston Rd
 Sharston
 Manchester M22 4SN
 Great Britain

Telephone 01382 826620
 Telefax 01382 826641
 E-mail: sales@juzo.com